



PPE International Trading

Merwedeweg 11
3336 LG Zwijndrecht
The Netherlands

+31 (0)85 029 0149
info@ppe-it.com
www.ppe-it.com

Model: ZSFM-01

FFP2 NR

CE 2834

USER INSTRUCTIONS


- 1 Failure to follow all instructions and limitations could seriously reduce the effectiveness of this particle half mask and could lead to illness, injury or death.
- 2 A properly selected respirator is essential, before occupational use, the wearer must be trained by the employer in the correct use of the respirator in accordance with applicable safety and health standards.
- 3 This particle half mask does not supply oxygen. Use it in adequately ventilated and where containing sufficient oxygen to support life.
- 4 Discard the respirator and replace with a new one if:
 - a. Excessive clogging of the particle half mask cause breathing difficulty.
 - b. The particle half mask become damaged
- 5 Leave the contaminated area if dizziness, irritation or other distress occurs

USE LIMITATIONS

- 1 Do not use the respirator to enter or stay in a contaminated area under the following circumstance:
 - a. Atmosphere contains less than 19.5% oxygen.
 - b. If you smell or taste contaminant
 - c. For protection against gases or vapors
 - d. Contaminants or their concentrations are unknown or immediately dangerous to life or health.
 - e. For sandblasting, paint-spray operations and asbestos.
 - f. In explosive atmospheres
- 2 Do not modify or misuse the mask
- 3 Do not use the particle half mask with facial hair or any other conditions that may prevent a good face-seal, the requirements for leakage will not be achieved
- 4 This respirator helps protect against certain particulate contaminants but does not eliminate exposure to out the risk of contracting disease or infection. Misuse may result to sickness or death
- 5 This particle half mask marked "NR", shall not be used for more than one shift

FITTING INSTRUCTION

- 1 Place fingers inside the particle filtering half mask. Bend the nose clip around fingers to form nose shape
- 2 Hold the particle filtering half mask in position over the nose and mouth, pull the headharness behind the ears.
- 3 Attach the headharness to the retaining clip, improve comfort and prevent leakage.
- 4 Use both hands, form the nose clip around the nose and shift the headharness to ensure a positive face seal.
- 5 Seal check
 - a. To test the fit of the particle filtering half mask without an exhalation valve, cup both hands over particle filtering half mask and inhale sharply.
 - b. If air flow is felt in the nose area, re-adjust and tighten the nose clip.
 - c. If air flow is felt around the edges of the respirator, adjust the particle filtering half mask /headharness to achieve a better fit
- 6 Change the particle filtering half mask immediately if breathing becomes difficult or particle filtering half mask becomes damaged or distorted
- 7 Change the particle filtering half mask if proper face seal cannot be achieved
- 8 Carefully observe these instructions is an important step in safe respirator use



Temperature range of storage conditions.



< 80%

Maximum relative humidity of storage conditions



See information supplied by the manufacturer.



yyyy/mm

Filtering half masks meet the requirements of Regulation (EU) 2016/425. They have been tested according to EN 149:2001+A1:2009 and certification by CCQS Certification Services Limited, Block 1 Blanchardstown Corporate Park, Ballycoolin Road, Blanchardstown, Dublin 15, D15 AKK1, Ireland (NB 2834).

Manufacturer: Wuhan Zonsen Medical Products Co., Ltd.

Address: Chengbei Industry Zone, Zhucheng Avenue, Xinzhou District, Wuhan City, China